


Raman microspectroscopy of MoS<sub>2</sub> by uRaman

In Yee Phang Ph.D, Eddie K. M. Tan, Ph.D, Technospex Pte. Ltd., Singapore

## Introduction

A monolayer of molybdenum disulfide (MoS<sub>2</sub>) has a unique and strong photoluminescence (PL) in the visible wavelength because of its direct bandgap at 1.8 eV (~680 nm). This unique characteristic makes monolayer MoS<sub>2</sub> a good candidate for optoelectronic application. This is unlike bulk MoS<sub>2</sub> with a few layers of MoS<sub>2</sub>, which is an indirect bandgap semiconductor with reduced PL.

Raman spectroscopy is a powerful technique to assess the number of layer of 2D materials such as in graphene and MoS<sub>2</sub>. In MoS<sub>2</sub>, the peak separation of in-plane S-Mo-S vibration ( $E_{2g}^1$  mode) and out-of-plane vibration of S atoms ( $A_{1g}$  mode) can be used to estimate the number of layers [1,2]. We measured the MoS<sub>2</sub> samples using uRaman from Technospex on a Nikon-Ci body. The uRaman module consists of a compact narrow line width laser (532 nm excitation and spectra bandwidth < 1 MHz). Using diamond as a reference, the system spectral resolution was measured to yield ~ 8 cm<sup>-1</sup>. Nikon 100X objective was used to measure the MoS<sub>2</sub> samples (model - TU Plan Fluor 100X, NA 0.9). The exfoliated MoS<sub>2</sub> sample was prepared by peeling a bulk MoS<sub>2</sub> using Scotch tape. By performing Raman spectroscopy on the exfoliated sample (Figure A), we have identified 6 different locations having different peak separations, indicating the presence of 6 different layers of MoS<sub>2</sub> (Figure B). The two spectra fingerprint,  $E_{2g}^1$  and  $A_{1g}$  bands, are observed at 383 and 401 cm<sup>-1</sup> at location 1. As the number of layers increases, these two peaks become more separated, reaching 380 and 405 cm<sup>-1</sup> at location 6 (Figure B and C). The Raman peak separation is contributed by the change in chemical interactions within the MoS<sub>2</sub> layers. As the layer number increases, the MoS<sub>2</sub> layers experience higher van der Waals force and long range Columbic interlayer interactions, which contributed to the blue shift in  $E_{2g}^1$  mode and red-shift in  $A_{1g}$  peak.



**Acknowledgement-** We thank Mr. Calvin Pei Yu Wong (NUS) for providing the MoS<sub>2</sub> sample for measurement.

Reference: [1] ACS Nano 2010, 4, 2695; [2] Adv. Funct. Mater., 2012, 22, 1385.