



NIUMAG

No Cryogen MRI  
Zero Maintenance



## Introducing the mini EDU-MR

The perfect solution for teaching MRI

### Features:

Easily accommodated on desktop  
or benchtop

### Specifications:


Footprint: 30W × 20D × 20H

Total Weight: 110 Jbs

Power: Standard

Sample size: 10mm

Field Strength: 0.5T



Time resolved images of water permeation  
into a white bean over a 24 hour period at  
room temperature

Suzhou Niumag Corporation

Web: [www.nmranalyzer.com](http://www.nmranalyzer.com) Tel: +86-512-62393560 Email: [info@niumag.com](mailto:info@niumag.com)

Korea South

Web: [www.jkxtac.com](http://www.jkxtac.com) Tel: +82-2-2145-4360 Email: [jkbhang@jktac.co.kr](mailto:jkbhang@jktac.co.kr)

# mini EDU-MR




Spatially Resolved Phantom

Contrast Media MRI

## Key Features:

- Durable for hands-on student use
- Virtual Software for student preparation
- Multiple imaging protocols
- Laptop or desktop PC to acquire data
- Easy to use software to tailor experiments



Sesame Seeds



Chile



Chicken Kunckle



Peanut in Shell